




Terumo® Oral / Enteral Syringes



Safe connection for oral and enteral administration

- Coloured plunger for visual differentiation
- Clear indication of intended use on barrel
- Special tip design for preventing wrong routing
- Long tip for easy oral administration
- NPSA conform



 **TERUMO®**

We keep life flowing

The special tip design of Terumo's oral and enteral syringes prevents misconnection to standard intravenous accessories as the diameter of the oral and enteral syringe's tip is too wide and consequently helps to avoid possible misconnections.



A coloured plunger for visual differentiation enables the easy recognition of the syringe type. Furthermore, it may equally allow you to integrate the use of oral and enteral syringes in existing procedures by allocating its colour to its intended use.

The clear indication of the intended use on the barrel gives your medical staff an additional visual confirmation and helps them to prevent a possible wrong routing.

The long oral tip can easily be inserted into the patient's mouth and thus minimizing the spilling of any medication.

Terumo® - Oral Syringe

Product Code	Description		Packaging	
	Size (ml)	Type	Unit box	Shipping carton
BS-50CO	50ml	Wide luer slip tip	25	100

Shelf life: 5 years



Sterilized by ethylene oxide

Terumo® - Enteral Syringe

Product Code	Description		Packaging	
	Size (ml)	Type	Unit box	Shipping carton
SS+50EO	50ml	Female luer lock	25	300
SS+30EO	30ml	Female luer lock	35	560

Shelf life: 5 years



Sterilized by ethylene oxide

Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients. Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



We keep life flowing

TERUMO CORPORATION
44-1, 2-chome
Hatagaya
Shibuya-ku
Tokyo 151-0072
JAPAN
Tel. +81/3 3374.8111
Fax +81/3 3374.8196

TERUMO EUROPE N.V.
Researchpark Zone 2
Haasrode
Interleuvenlaan 40
B-3001 Leuven
BELGIUM
Tel. +32/16.38.12.11
Fax +32/16.40.02.49

TERUMO EUROPE N.V. BENELUX SALES DIVISION
Researchpark Zone 2
Haasrode
Interleuvenlaan 40
B-3001 Leuven
BELGIUM
Tel. +32/16.39.25.80
Fax +32/16.39.25.99
THE NETHERLANDS
Tel. 0800/022.03.96
Fax 0800/022.04.14

LABORATOIRES TERUMO FRANCE S.A.
1, Parc Ariane
Bâtiment Uranus
Rue Hélène Boucher
78284 Guyancourt Cedex
FRANCE
Tel. +33/1 30.96.13.00
Fax +33/1 30.43.60.85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87
65760 Eschborn
GERMANY
Tel. +49/61.96.80.230
Fax +49/61.96.80.23.200

TERUMO DEUTSCHLAND GmbH
Zweigniederlassung
Wettingen
Bahnhofstrasse 88
CH-5430 Wettingen
SWITZERLAND
Tel. +41/56.419.10.10
Fax +41/56.419.10.11

TERUMO EUROPE ESPAÑA SL
Avda. Juan Carlos I,
N° 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares
(Madrid)
SPAIN
Tel. general:
+34/902.10.12.98
Tel. cust. serv.:
+34/902.20.31.40
Fax general:
+34/902.10.13.58
Fax cust. serv.:
+34/902.10.17.57

TERUMO ITALIA SRL
Via Simone Martini
143/145
I-00142 Roma
ITALY
Tel. +39/06.519.61.420
Fax +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
SE-426 71 Västra
Frölunda
SWEDEN
Tel. +46/31.748.58.80
Fax +46/31.748.58.90

www.terumo.com
www.terumo-europe.com

TERUMO DENMARK Filial af TERUMO SWEDEN AB
Tel. +45/7020.93.80
Fax +45/7020.94.80

TERUMO UK Ltd.
Tamesis
The Causeway
Egham
Surrey
TW20 9AW
UNITED KINGDOM
Tel. +44/1784.476.200
Fax +44/1784.476.201

TERUMO EUROPE N.V.
Moscow Representative office
4, Ivana Franko St
121108 Moscow
RUSSIA
Tel. +7/495.234.3303
Fax. +7/495.234.3303

®: Registered Trademark
Published by Terumo Europe N.V.

Addresses are subject to change, please consult our website.