





Terumo® Spinal Needle



Patient friendly penetration

- Smooth penetration
- Smooth cutting edge
- Unbendable and unbreakable needle
- Ergonomic, finger-fitting translucent grip-hub
- Easy orientation
- Total colour-coding



A smooth penetration with an optimal penetration resistance thanks to the spinal needle's bevel, which is oriented at a 22° angle. This angle gives you an optimal feel on where exactly the needle point is.

The smooth and rounded cutting edges of the cannula and stylet points are perfectly aligned. Due to the precise fit, there is minimal friction between both parts and the stylet can move easily in and out the cannula.

The needle will not bend or break, not even in case of an unexpected patient spasm as it is made of special high tensile materials. The danger of an inadvertent detachment is non existing, as both cannula and stylet are firmly affixed to the hub.

Designed to fit to your fingertips the grip-hub will ensure that your fingers are not cramped and you can push the needle without straining. The grip-hub will equally allow you to accurately confirm the flow of spinal fluid.

To easily determine the orientation reading of the needle point, a TERUMO logo is embossed, aligned with the needle bevel.

A total colour-coding is used on the stylet-hub, the blister pack and the shelf package for rapid gauge identification at every stage.



Terumo® Spinal needle

Product Code	Colour Code		Dimensions of cannula		
			Diameter		Length
SN*1870		Pink	18G	(1.2 mm)	70 mm
SN*1890		Pink	18G	(1.2 mm)	90 mm
SN*1970		Brown	19G	(1.1 mm)	70 mm
SN*1990		Brown	19G	(1.1 mm)	90 mm
SN*2038		Yellow	20G	(0.9 mm)	38 mm
SN*2070		Yellow	20G	(0.9 mm)	70 mm
SN*2090		Yellow	20G	(0.9 mm)	90 mm
SN*2138		Green	21G	(0.8 mm)	38 mm
SN*2170		Green	21G	(0.8 mm)	70 mm
SN*2190		Green	21G	(0.8 mm)	90 mm
SN*2238		Black	22G	(0.7 mm)	38 mm
SN*2270		Black	22G	(0.7 mm)	70 mm
SN*2290		Black	22G	(0.7 mm)	90 mm
SN*2370		Light blue	23G	(0.6 mm)	70 mm
SN*2390		Light blue	23G	(0.6 mm)	90 mm
SN*2590		Blue	25G	(0.5 mm)	90 mm

Sterilized by gamma.
Unit box: 20 pieces
Shipping carton: 200 pieces
Sterilized by ethylene oxide



Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients. Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



TERUMO CORPORATION
44-1, 2-chome
Hatagaya
Shibuya-ku
Tokyo 151-0072
JAPAN
Tel. +81/3 3374.8111
Fax +81/3 3374.8196

TERUMO EUROPE N.V.
Researchpark Zone 2
Haasrode
Interleuvenlaan 40
B-3001 Leuven
BELGIUM
Tel. +32/16.38.12.11
Fax +32/16.40.02.49

TERUMO EUROPE N.V. BENELUX SALES DIVISION
Researchpark Zone 2
Haasrode
Interleuvenlaan 40
B-3001 Leuven
BELGIUM
Tel. +32/16.39.25.80
Fax +32/16.39.25.99
THE NETHERLANDS
Tel. 0800/022.03.96
Fax 0800/022.04.14

LABORATOIRES TERUMO FRANCE S.A.
1, Parc Ariane
Bâtiment Uranus
Rue Hélène Boucher
78284 Guyancourt Cedex
FRANCE
Tel. +33/1 30.96.13.00
Fax +33/1 30.43.60.85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87
65760 Eschborn
GERMANY
Tel. +49/61.96.80.230
Fax +49/61.96.80.23.200

TERUMO DEUTSCHLAND GmbH Zweigniederlassung
Wettingen
Bahnhofstrasse 88
CH-5430 Wettingen
SWITZERLAND
Tel. +41/56.419.10.10
Fax +41/56.419.10.11

TERUMO EUROPE ESPAÑA SL
Avda. Juan Carlos I,
N° 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares
(Madrid)
SPAIN
Tel. general:
+34/902.10.12.98
Tel. cust. serv.:
+34/902.20.31.40
Fax general:
+34/902.10.13.58
Fax cust. serv.:
+34/902.10.17.57

TERUMO ITALIA SRL
Via Simone Martini
143/145
I-00142 Roma
ITALY
Tel. +39/06.519.61.420
Fax +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
SE-426 71 Västra
Frölunda
SWEDEN
Tel. +46/31.748.58.80
Fax +46/31.748.58.90

TERUMO DENMARK Filial af TERUMO SWEDEN AB
Tel. +45/7020.93.80
Fax +45/7020.94.80

TERUMO UK Ltd.
Tamesis
The Causeway
Egham
Surrey
TW20 9AW
UNITED KINGDOM
Tel. +44/1784.476.200
Fax +44/1784.476.201

TERUMO EUROPE N.V. Moscow Representative office
4, Ivana Franko St
121108 Moscow
RUSSIA
Tel. +7/495.234.3303
Fax. +7/495.234.3303

®: Registered Trademark
Published by Terumo Europe N.V.



We keep life flowing

Addresses are subject to change, please consult our website.

www.terumo.com
www.terumo-europe.com