

Surflo® Winged Infusion Set



Perform and Care

- *Smooth placement*
- *Ease of use*
- *Patient comfort*
- *User satisfaction*

Surflo® Winged Infusion Set

Smooth placement

Terumo® cannula technology ensures an outstanding sharpness combined with a special silicone coating allowing an easy and smooth placement.

Ease of use

The user inspired design of the large and interlocking wings ensures a firm grip during placement. Colour coded wings and kink resistant tubings provide additional ease of use.


Patient comfort

The excellent quality of the Terumo® cannula is at the core of minimising trauma. Patients reported having experienced an almost painless placement.


User satisfaction

The design and performance of the Surflo® winged infusion sets are supporting a successful placement and increasing end-user satisfaction due to less discomfort caused to patients.


Instructions for use



After successful infusion, hold the winged needle in place.



Withdraw the winged infusion set carefully from the vein.



Dispose the winged infusion set after use into a disposal box.


Products and specifications

Code	Colour	Cannula - external Ø		Cannula length	Tubing length	Dead space mL	Packaging
		G	mm				
SV-18BL	Pink	18	1.2	19 mm	300 mm	0.44	10 x 50
SV-19BL	Ivory	19	1.1	19 mm	300 mm	0.43	10 x 50
SV-21BL	Green	21	0.8	19 mm	300 mm	0.42	10 x 50
SV-22BL	Black	22	0.7	19 mm	300 mm	0.42	10 x 50
SV-23BL	Blue	23	0.6	19 mm	300 mm	0.41	10 x 50
SV-25BL	Orange	25	0.5	19 mm	300 mm	0.41	10 x 50
SV*27EL	Grey	27	0.4	13 mm	200 mm	0.27	10 x 100

Shelf life: 5 years; 3 years (SV*27EL)



CE



Regular Wall

Thin Wall

Terumo-Thin Wall



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

European Offices

**Terumo Europe N.V.
Benelux Sales Division**
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.
Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex, France
Tel.: +33 130 96 13 00
Fax: +33 130 43 60 85

Terumo Deutschland GmbH

Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Europe España SL
Avda. Juan Carlos I, Nº 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares (Madrid),
Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Italia S.r.l.
Via Simone Martini, 143/145
00142 Roma RM,
Italy
Tel.: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo Sweden AB

Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
426 71 Västra Frölunda,
Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

**Terumo Denmark
Filial of Terumo Sweden AB**
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

Terumo UK Ltd
Tamesis
The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 200
Fax: +44 1784 476 201

**Terumo Europe N.V.
Moscow Representative Office**
BC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: +7 495 988 4740
Fax: +7 495 988 4739

*Addresses subject to change,
please consult our Website.*

© Registered Trademark
Published by Terumo Europe N.V.
HO 07GB-0406FK-III(08.11)E