




Terumo[®] Syringes

Two part syringes



Controlled aspiration and injection for medical staff

- Highly transparent barrel for visualisation
- Firm and definite end stop
- Ergonomic finger grips
- Latex and PVC-free
- Clear scale markings to enable precise dosage



 **TERUMO[®]**

We keep life flowing

The ergonomically designed finger grips and plunger push-button allow a perfectly controlled and stable injection.

Firm and definite end-stop prevents accidental separation or spillage, avoiding exposure to hazardous fluids.

The latex and PVC-free syringes are made of carefully selected, environmental friendly materials and are being systematically controlled during the production process.

The highly transparent barrel allows for easy recognition of air bubbles and the clear scale markings allow a precise dosage management.

The sterile, individual peel-off blister pack of each syringe keeps it in a sterile condition until it will be used. The lot number and expiry date are clearly marked on each blister.



Terumo® two part syringes

Product Code	Description		Packaging	
	Size (ml)	Type	Unit box	Shipping carton
SS+T02S1	2.5	Centered tip	100	2800
SS+T05ES1	5	Eccentric tip	100	1900
SS+T10ES1	10	Eccentric tip	100	1200
SS+T20ES1	20	Eccentric tip	50	700

Sterilized by ethylene oxide



Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients. Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



TERUMO CORPORATION
44-1, 2-chome Hatagaya Shibuya-ku Tokyo 151-0072 JAPAN
Tel. +81/3 3374.8111 Fax +81/3 3374.8196

TERUMO EUROPE N.V.
Researchpark Zone 2 Haasrode Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.38.12.11 Fax +32/16.40.02.49

TERUMO EUROPE N.V. BENELUX SALES DIVISION
Researchpark Zone 2 Haasrode Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.39.25.80 Fax +32/16.39.25.99
THE NETHERLANDS
Tel. 0800/022.03.96 Fax 0800/022.04.14

LABORATOIRES TERUMO FRANCE S.A.
1, Parc Ariane Bâtiment Uranus Rue Hélène Boucher 78284 Guyancourt Cedex FRANCE
Tel. +33/1 30.96.13.00 Fax +33/1 30.43.60.85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87 65760 Eschborn GERMANY
Tel. +49/61.96.80.230 Fax +49/61.96.80.23.200

TERUMO DEUTSCHLAND GmbH Zweigniederlassung Wettingen
Bahnhofstrasse 88 CH-5430 Wettingen SWITZERLAND
Tel. +41/56.419.10.10 Fax +41/56.419.10.11

TERUMO EUROPE ESPAÑA SL
Avda. Juan Carlos I, N° 13 - 7ª Planta Edificio Torre La Garena 28806 Alcalá de Henares (Madrid) SPAIN
Tel. general: +34/902.10.12.98
Tel. cust. serv.: +34/902.20.31.40
Fax general: +34/902.10.13.58
Fax cust. serv.: +34/902.10.17.57

TERUMO ITALIA SRL
Via Simone Martini 143/145 I-00142 Roma ITALY
Tel. +39/06.519.61.420 Fax +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90 Sven Källfelts Gata 18 SE-426 71 Västra Frölunda SWEDEN
Tel. +46/31.748.58.80 Fax +46/31.748.58.90

TERUMO DENMARK Filial af TERUMO SWEDEN AB
Tel. +45/7020.93.80 Fax +45/7020.94.80

TERUMO UK Ltd.
Tamesis The Causeway Egham Surrey TW20 9AW UNITED KINGDOM
Tel. +44/1784.476.200 Fax +44/1784.476.201

TERUMO EUROPE N.V. Moscow Representative office
4, Ivana Franko St 121108 Moscow RUSSIA
Tel. +7/495.234.3303 Fax. +7/495.234.3303

®: Registered Trademark Published by Terumo Europe N.V.



Addresses are subject to change, please consult our website.

www.terumo.com
www.terumo-europe.com